



### Low-temperature liquid helium sample holder for cryo-TEM

- Variable temperature
- Continuous helium flow
- Cool-down 10– 15 min
- 16 electrical contacts

### Milestones

- Founded in 2019
- Patent filed for novel cryostat technology
- UZH Bio/MedTech Entrepreneur Fellowship

### Expertise

Development and production of scientific instruments and components for research applications in ultra-high vacuum and cryogenic conditions.


### Contact

condenZero GmbH  
 % Physik-Institut  
 Universität Zürich  
 Winterthurerstrasse 190  
 8057 Zürich  
 +41 44 635 57 80


### Team



Dr. Denys Sutter



Dominik Biscette



Prof. Johan Chang